



eQuinox™

Transportable Oxygen Concentrator with autoSAT®

The perfect balance between portability and clinical efficiency

The SeQual eQuinox is the newest, quietest, and most advanced transportable POC on the market. With cutting-edge technology providing numerous new features to meet both patient wants and provider needs, the SeQual eQuinox is the premier POC available.

The eQuinox continues to deliver the revolutionary clinical benefits you have come to expect from SeQual POCs, but now includes additional user enhancements. The Multi-Language Voice Interface offers a new layer of comfort to users by providing verbal confirmation of changed flow rate settings, as well as announcing battery times and any alarms. The language can even be changed to offer feedback in the language most comfortable to the end user.

Weighing in at only 14 pounds, the newly designed, easy-to-manuever frame makes the SeQual eQuinox the lightest weight POC able to offer 3.0 LPM continuous flow. This small package is powerful enough to give patients the freedom to enjoy life on the go. The eQuinox is fully functional on AC Power, DC Power, and battery power. Wherever a patient goes, the eQuinox goes with them and keeps up with the various activities of daily living. With the option of a longer lasting 24-cell battery pack, the eQuinox can achieve 2.75 hours of operation at 2.0 LPM.

The eQuinox is an all-in-one non-delivery modality solution for your oxygen business. Providing continuous flow options from 0.5LPM to 3.0LPM and 9 pulse flow settings from 16mL to 192mL, the SeQual eQuinox can meet a patient's 24-7 needs. The comprehensive dosing selection allows a patient to be treated at rest, during sleep, at exercise, and at altitude.

Features

- **NEW: Multi-Language Voice Interface makes using a POC simpler than ever**
- **NEW: Only 14 lbs with 12-cell battery pack installed**
- **NEW: Up to 2.75 hours of battery life at 2.0 LPM with the optional 24-cell battery pack**
- **NEW: Only 45 dB at 2 continuous flow**
- **NEW: Service Cabinet makes filters and the 9V battery accessible for service technicians without removing the case—performing scheduled maintenance and upkeep is easier than ever**
- **NEW: DC Power supply consists of a lightweight 12V cable—no transformer brick**
- **NEW: Screen intermittently displays battery charge remaining in minutes**
- Approved for RX of 0.5 - 3.0 LPM Continuous Flow
- 9 pulse settings with bolus sizes from 16mL to 192mL
- autoSAT technology helps maintain a consistent FIO2 by servo-controlling the device to meet changing respiratory rates
- Three available rise time options: Fast, Medium, Slow
- Three adjustable Pulse Dose Sensitivity settings with a maximum sensitivity of -0.13 cm H2O
- FAA Approved for use on commercial aircraft
- Full functionality on DC power
- Ability to recharge battery cartridge on DC power when excess power is available



Call 0870 712 0202 to order!

www.healthoxygen.com

eQuinox™

Transportable Oxygen Concentrator with autoSAT®

Product Specifications



	eQuinox Oxygen System
Size H x W x D	13.6 x 10.6 x 7.4 (inches)
Weight	14 lbs. (6.4 kg) with 12-Cell Battery Pack Installed
Continuous Flow Settings	0.5 to 3.0 LPM (0.5 liter increments)
Pulse Dose Settings	16-96mL (16mL increments), 128mL, 160mL, 192mL
O ₂ concentration	90%+5.5%/-3% @ sea level
O ₂ output pressure (nominal)	5.0 psig (34.5 kPa)
Sound level	46dBA @ 3.0 LPM Continuous Flow; 45dBA @ 2.0 LPM Continuous Flow; 40dBA @ 0.5 LPM Continuous Flow
Power AC Power	85-264VAC, 47-63Hz
DC Power	12 V nominal
Battery Pack	Quantity (1) 89 W-Hr battery pack 12-cell battery; Quantity (2) 95 W-Hr battery packs 24-cell battery
Power Consumption	110 watts at 3.0 LPM Continuous
Battery Duration	1.2 hrs at 2 LPM Continuous; 2.7 hrs at 32mL pulse (Setting 2) (12 BPM) with 12-cell Battery Pack 2.75 hrs at 2 LPM Continuous; 5.9 hrs at 32mL pulse (Setting 2) (12 BPM) with optional 24-cell Battery Pack
Nominal Battery Recharge Time	1.5 to 3.5 hrs recharge time to achieve 80% capacity (dependent upon the flow)
Backup Alarm power	9 volt battery
O ₂ concentration indicator	Green Light=Normal, Yellow=Caution <85%, Red Flashing=Warning <70%
Operating temperature	32°F to 104°F (0°C to 40°C) at 82.4°F (28°C) dew point, non-condensing
Trigger Sensitivity	Adjustable between -0.13 cm H ₂ O to -0.4 cm H ₂ O (negative pressure)
Altitude Operating range	-382 to 4000 m International Standard Atmosphere (-1254 ft to 13,130 ft) (1060 to 616 hPa)
Alarms/Alerts	Loss of power, Low Battery, Low Therapeutic O ₂ Output, O ₂ No Inspiration detected in Pulse Dose Mode, Unit Malfunction

AIRLINE APPROVALS

eQuinox is approved by the Federal Aviation Administration (FAA) for onboard, in-flight use.

Ordering Information

Description	Part Number
eQuinox Oxygen System consists of Concentrator, Battery Pack, AC Power Supply, DC Power Cable, Accessory Bag and Cart with telescoping handle.	4807-SEQ

Pure O₂

Viking House, 71 Princess Road, Urmston, Manchester M41 5ST

Call 0870 712 02 02

Website: www.healthoxygen.com Email: info@healthoxygen.com

Pure O₂ reserves the right to discontinue its products, or change the prices, materials, equipment, quality, descriptions, specifications and/or, processes to its products at any time without prior notice and with no further obligation or consequence. All rights not expressly stated herein are reserved by us, as applicable.

Copyright © 2014 Chart Industries ML-CONC0033 Rev. A