

# Pure O2 introduce the NEW, INOGEN ONE® G3 5lpm

**SIMPLY THE BEST PORTABLE OXYGEN CONCENTRATOR AVAILABLE**

*Pulse settings up to 5lpm of oxygen delivery whilst weighing only 2.2kg and running very quiet at 39 Decibels*



Only  
4.8  
pounds!

## The Inogen One® G3 System (IS-300):

- Compact and Lightweight
- Noticeably Quiet Operation
- Intelligent Oxygen Delivery
- Designed for 24/7 Use

\*Included with system

\*\*Accessory sold separately

- An additional 5th setting - now serve an even wider patient group;
- 1050 ml/min of oxygen delivery - highest output on the market;
- Sound level of only 39 dBA - quietest POC compared to any competitive brand;
- New gross particle filters - more robust and easier to clean.

pureO<sub>2</sub>



Inogen One®  
AC Power Supply\*  
Item # BA-301



Inogen One® G3  
DC Power Cable\*  
Item # BA-306



Inogen One® G3  
Carry Bag\*  
Item # CA-300



Inogen One® G3  
Lithium Ion Battery\*  
Up to 4.5 hours run time  
Item #BA-300



Inogen One® G3  
Extended Life Lithium Ion Battery\*\*  
Up to 9 hours run time  
Item #BA-316



INOGEN ONE® G3

## OXYGEN. ANYTIME. ANYWHERE®

The Inogen One G3 delivers the independence of a portable oxygen concentrator in one of the smallest, lightest, and quietest packages available to the oxygen user today. With flow settings from 1-5, the Inogen One G3 is designed to increase independence for most users 24/7. It can be charged at home or on the go - giving you the freedom of Oxygen. Anytime. Anywhere.®

The Inogen One® G3 is your one solution, your single solution, for oxygen at home or away.

<b>Weight</b>	4.8 lbs. (includes single battery)
<b>Size</b>	<b>Length:</b> 8.75" <b>Width:</b> 3.0" <b>Height:</b> 8.25" (includes single battery)
<b>Oxygen Flow</b>	Pulse dose delivery system Five flow settings 1, 2, 3, 4, 5
<b>Power</b>	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car
<b>Battery</b>	Duration (single battery): Up to 4.5 hours Recharge approximately 2-4 hours with AC or DC power Duration (double battery): Up to 9 hours Recharge approximately 4-8 hours with AC or DC power
<b>Noise Level</b>	39 Decibels*
<b>Warranty</b>	3 year or lifetime limited warranty**
<b>Operation</b>	Simple control functions and easy-to-read LCD display
<b>Use</b>	Designed for 24/7 use

\*At flow setting 2

\*\*See Inogen One® Limited Warranty Statement for complete warranty information

Call us now on **0161 747 2617**, our 24 hour line **0870 712 0202**  
or email us on: [info@healthoxygen.com](mailto:info@healthoxygen.com)

Supplied by **Pure O2 Ltd**  
Viking House, 71 Princess Road  
Urmston, Manchester  
M41 5ST, United Kingdom

[www.healthoxygen.com](http://www.healthoxygen.com)

