

# FREESTYLE<sup>®</sup> COMFORT<sup>™</sup>

Portable Oxygen Concentrator

## The future is here

With its ergonomic shape, the lightweight and affordable FreeStyle<sup>®</sup> Comfort<sup>™</sup> rests against the natural curves of your body, making it comfortable to wear on the go. FAA approved for commercial air flights and equipped with user-replaceable batteries, this new oxygen therapy device allows you to extend your time away from home. And extended warranty options offer peace of mind. It's time to expect more from your portable oxygen concentrator.



### A LIGHTWEIGHT SOLUTION

Weighing only 5 lbs., the FreeStyle<sup>®</sup> Comfort<sup>™</sup> can go anywhere you want to go.


Enjoy the convenience of wearing it with or without your carrybag by simply attaching the no-slip, padded shoulder strap.



### BREATHE EASY, WORRY LESS

UltraSense<sup>™</sup>, our sensitive breath detection technology, ensures oxygen is delivered to you when you need it.

User-replaceable batteries allow you to extend your time away from home and stay active with up to 4 hours of battery life.



### COMFORTABLE TO WEAR, EASY TO USE

The ergonomic, curved silhouette rests against the natural curves of your body making it more comfortable to wear when you are out and about.

An easy-to-use full-color interface will raise the bar on what you expect from your portable oxygen concentrator.



Contact Information:  
Distributor Pure O2  
Tel UK 0161 747 2617  
Int. 0044 161 747 2617  
[www.healthoxygen.com](http://www.healthoxygen.com)  
[www.cairemedical.com](http://www.cairemedical.com)

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.  
Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2018 CAIRE Inc. All Rights Reserved.

# FREESTYLE<sup>®</sup> COMFORT<sup>™</sup>

Portable Oxygen Concentrator

## Includes

- Portable Concentrator
- Carry Bag
- Battery Pack
- AC/DC Universal Power Supply and Cords
- Shoulder Strap
- UltraSense<sup>™</sup> technology
- Other SKU's may include extra accessories

## Product Features

- User-replaceable battery pack with up to 4 hours of battery life at setting 2
- Easy to wear with a shoulder strap
- FAA-approved for commercial airline use
- Fully functional on wall (AC) or motor vehicle (DC) power
- Easy-to-use full-color interface

## Specifications

### FLOW RATE

Pulse 1-5

### DIMENSIONS

10.0 x 7.3 x 3.1 in  
(25.4 x 18.5 x 7.9 cm)

8-cell battery pack

11.0 x 7.3 x 3.1 in

(27.9 x 18.5 x 7.9 cm)

16-cell battery pack

### WEIGHT

5 lbs (2.3 kg)

8-cell battery pack

6 lbs (2.7 kg)

16-cell battery pack

### SOUND LEVEL

43 dB(A) @ Setting 2

### BATTERY DURATION

8-cell battery at setting  
2 – up to 4 hours

16-cell battery at setting  
2 – up to 8 hours

### WARRANTY

3 years FreeStyle<sup>®</sup> Comfort<sup>™</sup>  
(Includes sieve beds)

1 year battery and power supplies



## Available Accessories



**8-CELL  
BATTERY PACK**



**16-CELL  
BATTERY PACK**



**CARRYING BAG**



**DESKTOP BATTERY  
PACK CHARGER  
(OPTIONAL)**



**AC POWER SUPPLY  
(ALL CORDS INCLUDED)**



**DC POWER CORD**



**EXTERNAL GROSS  
PARTICLE FILTERS**

### OTHER ITEMS AVAILABLE:

- Shoulder Strap
- Carry-all Accessory Bag



Distributor Pure O2  
Tel UK 0161 747 2617  
Int. 0044 161 747 2617  
[www.healthoxygen.com](http://www.healthoxygen.com)

[www.cairemedical.com](http://www.cairemedical.com)

SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.  
Please consult the applicable product instructions for use for product indications, contraindications, warnings, precautions, and detailed safety information.

© 2018 CAIRE Inc. All Rights Reserved.

ML-CONC0138 C